

ATF Video Frame Grabber Subsystems

- Frequently Asked Questions
- Tips and Techniques

For ATF Staff, Operators and Visiting Users

Contact: R. Malone, X2084; malone@bnl.gov

Last updated: Tue Aug 4 15:21:48 EDT 2009
Revision history at end of document.

=====
Part II - Saving and recalling images and their associated color maps

(Q1) Where are the image and color map files stored?

(A1) By default, on the frame grabber PC where you originally captured the data, in the "C:\Images" folder. Most users create a sub-directory under this top level folder, giving it a name related to their experiment and the date.

You are not restricted to using "C:\Images". You may use any device/folder for which you have write permission, including external devices such as your USB stick.

Image files are stored with a file type of ".asc" and color maps are stored with a file type ".atf_cm".

(Q2) How can I copy my saved image data from the frame grabber PC on to my own PC/USB stick/other media?

(A2) You can connect your USB device to the front panel USB connector on the frame grabber PC chassis. Once your device is mounted, you can click/drag files from the frame grabber's image folder to your external device.

OR

You can connect your PC to the local private ATF network. Set your PC to use DHCP and an IP address will be assigned to you. You can then start an ftp client on your PC, instructing it to connect to the frame grabber PC where you captured your images. Use the network node names as explained in Section I of this FAQ. This method may be more convenient if you are transferring a large amount of data.

Ask ATF staff for the username/password needed to access the systems via ftp.

(Q3) How can I reload a stored image file?
How can I verify that a stored image file is OK?
How can I reload a color map?
How can I verify that a stored color map is OK?

(A3) You can recall a stored image file into the ATF frame grabber application by pressing the "Recall" button and navigating to the device/folder where you stored the image. Select the image file and it will be read into the frame grabber application and will be treated as a newly acquired image. If the image displays correctly, you can be confident the image file was saved correctly.

Similarly, you can reload a color map from a stored file. The color bar at the right of the image area should change to reflect the newly loaded color map.

(Q4) What's stored in a captured image file and how is it arranged?

(A4) The captured image file is an ASCII text file which contains the raw, unprocessed video data as digitized by the frame grabber or GigeVision camera.

The data is stored in row-major order, proceeding from the top left of the image and moving vertically line by line, like this:

```
data for video line 1 ...
data for video line 2 ...
data for video line 3 ...
.
.
.
data for last video line.
```

Each line consists of comma separated values (CSV), with each value being the pixel intensity level. The last value on each line does not have a trailing comma. Each line is terminated with two ASCII characters: <CR> and <LF>.

All data captured by analog frame grabbers has been fully de-interlaced.

Sample fragment of image file:

```
6,11,13,3,3,9,9,10,7,11,1,3,5,2,4,5,8,0,11,14,6,3 ... <-- row 1 of image
10,7,6,6,0,19,3,1,0,11,2,7,0,19,16,4,1,4,5,5,7,12 ... <-- row 2 of image
0,1,1,4,7,8,2,6,7,3,9,0,8,0,14,0,8,14,5,17,5,11,0 ... <-- row 3 of image
.
.
.
11,11,6,8,7,0,3,0,3,0,1,0,20,4,0,2,2,6,5,2,3,10,9 ... <-- last row of image
```

(Q5) Why is only the raw, unprocessed data stored in an ATF image data file?

(A5) Having the original, raw data available for off-line analysis is valuable since it can be re-analyzed using different color maps, gain/black level, etc. If we stored the processed data, you would not be able to go back to the original data and re-process it.

(Q6) What are the dimensions of standard ATF image files?

(A6) The dimensions depend on both the type of frame grabber and camera in use when the image was captured:

For Frame Grabbers 5, 6, 7 and 8:

```
Camera type : Analog RS-170
Dimensions  : 480 rows x 640 columns
```

For Frame Grabbers 10, 11, 12:

```
Camera type : Basler Scout scA1400 17gm GigE
Dimensions  : 1040 rows x 1392 columns
```

Note: At present, the only GigE Vision cameras at ATF are the Basler Scouts with 1040 x 1392 image data. Future GigE Vision cameras may be added to ATF's inventory which most likely will have different dimensions. Any analytical software you plan to write and bring to ATF should be prepared to dynamically

determine the appropriate matrix size from the CSV image files.

- (Q7) Why do the dimensions I get when reading in an ATF image file differ from those specified by the camera manufacturer on their data sheets?
- (A7) Matlab (and other matrix/mathematically-oriented packages) declare matrix dimensions in the standard "rows x columns" form while camera manufacturers state image dimensions in "y-direction by x-direction" (i.e., "columns x rows") form.
-

(Q8) How is the color map stored?

- (A8) The color map is stored as an ASCII text file, with one line per possible pixel intensity level.

Each line consists of 3 comma-separated values in the closed range [0,1]. The last value in the line does not have a trailing comma. The line is terminated with two ASCII characters: <CR> and <LF>.

Each value represents the fraction of red, green and blue (RGB) needed to form the desired color.

Example: Sample fragment of 12-bit color map file

```
0,0,0.50098 <-- pixel intensity 0: Red 0, Green 0, Blue 0.50098
0,0,0.50195 <-- pixel intensity 1: Red 0, Green 0, Blue 0.50195
0,0,0.50293 <-- pixel intensity 2: Red 0, Green 0, Blue 0.50293
.
.
.
0.50195,0,0 <-- pixel intensity 4093: Red 0.50195, Green 0, Blue 0
0.50098,0,0 <-- pixel intensity 4094: Red 0.50098, Green 0, Blue 0
0.5,0,0 <-- pixel intensity 4095: Red 0.5, Green 0, Blue 0
```

(Q9) What are the dimensions of ATF color map files?

- (A9) The dimension of the data contained in an ATF color map file depend on the level of digitization which was in effect at the time the color map was captured:

```
For 8-bit digitization: 256 rows x 3 columns
For 10-bit digitization: 1024 rows x 3 columns
For 12-bit digitization: 4096 rows x 3 columns
```

(Q10) Why does my computer report different file sizes for image files I know I captured from the same frame grabber and at the same digitization level?

- (A10) Data items recorded in image capture files and color map files are not written in fixed-width formats. Instead, each data item is written using just enough characters to record the desired value without any leading or trailing blank/zero characters. The goal being, of course, to keep the file size as small as possible.

As an example, when digitizing with 12-bits, a 1040 x 1392 image captured by a Basler camera containing all zero pixels would require 2,894,320 bytes of storage, while one with fully saturated pixels would require 5,787,600 bytes.

(Q11) Is it important to capture and save the color map?

(A11) Maybe. If you find yourself needing to edit the default startup ATF color map so that you can bring out features in your image which would not otherwise be visible, you should probably save your color map after each significant edit. Keep a log of what color map file had been used with each captured image.

(Q12) What is the default ATF color map?

(A12) The default ATF color map is the Matlab color map named "jet", with the number of levels determined by the digitization level:

For 8-bit digitization, "jet(256)"
For 10-bit digitization, "jet(1024)"
For 12-bit digitization, "jet(4096)"

(Q13) In addition to capturing the raw image data to a file, is there any way to capture/print the current screen image so I can keep a copy with my lab notes?

(A13) Yes. Press the ALT and PRINT_SCREEN keys at the same time to bring up the screen capture utility (Print Screen 2000). Save the screen image in the desired format (.BMP, .JPG, etc.)

CAUTION: A screen capture can not take the place of a true image data file, particularly for any image analysis application. This is because pixels in the screen capture do not have a one-to-one correspondence with the pixels of the frame grabber image.