

EPICS Stream Device Programming

Marty Smith

APS Engineering Support Division

Agenda

- Introduction to Stream Device
- EPICS Databases and Stream Device
- Stream Device Protocols and Protocol Files
- Adding Stream Device Support to an Application
- Short Break
- Lab Session: Controlling a Network Attached Device



Introduction to Stream Device

- Generic EPICS device support for devices with a "byte stream" based communication
 - RS-232 (Local serial port or LAN/Serial adapter)
 - TCP/IP
 - VXI-11
 - GPIB (Local interface or LAN/GPIB adapter)
 - USB-TMC (Test and Measurement Class)
- A single stream device module can serve to communicate using any of the above communication mechanisms.

Introduction to Stream Device

- Command / Reply messages
 - *IDN?
 - xx:SetVoltageOut 1.2
 - Can include non-printable characters as well
- Command and reply parsing configured by **protocols**
- Formatting and interpretation handled with format converters
 - Similar to C printf and scanf format converters
 - Write your own converters too using the API



Introduction to Stream Device

- Provides interface to ASYN
 - Not a replacement for ASYN
 - ASYN driver provides string exchange to/from device
 - Uses AsynOctet interface
- Stream Device is not:
 - Programming language
 - No looping or conditionals available
 - Protocols are linear running from start to end
 - Rudimentary exception handlers
- How do we get stream device in a EPICS database?

Simple Command and Reply Message

- Simple command generating long response message

Data Sent: Q

Receive Data:

```
:SN=AT267 ,UN=id13 ,IP=164.054.008.127,V3=3390,V5=5135,V+12=12160,V-  
12=12396,T1=30,T2=28,T3=37,T4=00,F1=02160,F2=02130,F3=02160,F4=00000,F5=00  
000,F6=00000,F7=00000,F8=00000,F9=00000,OT=0,OV=0,OC=0000,PS1=1,PS2=1,MSG  
=0,SW=1,PROT=TEL ,I3=00,I5=00,I12=00,I-12=00,CODE=64-  
113426F39,ENET=D6.10,POH=28896.5,MAXTMP=43,MINTMP=22,PROC=31,LOAD=  
FF,PWRCYCL=00083  
;EV00000000ET00000000EF0000000000OT0OV0OC0000PS11MSG0SW1
```

- Protocol entries can be long

Stream Device EPICS Database

```
record(bo, "$($P)$($R)query") {  
    field(DESC, "Timed readback")  
    field(SCAN, "10 second")  
    field(PINI, "YES")  
    field(FLNK, "$($P)$($R)VP3")  
}  
record(ai, "$($P)$($R)VP3"){  
    field(DTYP, "stream")  
    field(INP, "@devDawnRuSH.proto query($($P)$($R)) $($PORT) 0")  
    field(FLNK, "$($P)$($R)VP5")  
}
```

- DTYP=stream
- INP/OUT fields specify protocol file name, protocol entry (with optional arguments), ASYN port and address.
- \$(\$P)\$(\$R) is a protocol argument, up to nine arguments can be provided
- Address can be any value (typically 0) for single-address interfaces

Stream Device EPICS Database

- DTYP \neq stream for protocol entry additional records:

```
record(stringin, "$(P)$ (R)Serial"){  
    field(DESC, "Serial number")  
    field(DTYP, "Soft Channel")  
}  
record(ai, "$(P)$ (R)VP5"){  
    field(DESC, "+5V supply")  
    field(DTYP, "Raw Soft Channel")  
    field(EGU, "V")  
    field(PREC, "3")  
    . . .
```


Stream Device Protocol Files

- Example protocol file
 - Use multiple lines to format for easy reading

```
query {  
    out "Q";  
    in ":SN=%(\$1Serial.VAL) 39[^,], "  
        "UN=%(\$1Name.VAL) 39[^,], "  
        "IP=%*[^,], "  
        "V3=%d, "  
        ...  
        "PWRCYCL=%(\$1PowerCycle.VAL)d";  
  
    ExtralInput = Ignore;  
}
```

- Notice the use of the width field – guard against buffer overruns!

Stream Device Protocols

- Defined in a plain ASCII text protocol file
- No need to compile, protocol is read by IOC at boot time
- A single entry can read/write multiple fields in one or many records
- Output records can be initialized from instrument at IOC startup
 - Providing instrument is powered on and communicating at IOC boot time
- Each protocol file describes communication to ONE device
- Protocols are defined for each function of a device

Stream Device Protocol Files

- All lines beginning with a # to the end of line are comments
- A protocol consists of a name followed by a body in {}
- Protocol entries contain statements to produce output and request input
 - Look similar to C functions
 - GetOutput {

```
        out “\ $1”;  
        in “%f”;  
    }
```
 - \$1 is a protocol argument, up to nine arguments can be provided

Stream Device Protocol Files

- By default the VAL or RVAL field is used as the data source/destination
- Can refer to any field, even in another record
- C-style escape sequence can be used ('\r', '\n', '\033', '\e')
- Can reload a protocol or all protocols without rebooting
 - Good for development of frequently changing protocol files
 - streamReload("recordname") – Reloads protocol for recordname
 - streamReload() – Reloads all protocols in a file

Stream Device Protocol Files

- Protocol file terminators
 - Terminators can be set globally or per entry
- Some interfaces can handle only a single character. If device replies with '\r\n' then specify InTerminator='\n' and ignore the '\r' in the reply
 - InTerminator = "\n";
 - OutTerminator = "\r";
- Better practice to use the ASYN terminators in IOC boot file
 - asynOctetSetOutputEOS and asynOctetSetInputEOS

Stream Device Protocol Files

- Initial read back from device at IOC boot time
 - Useful to set initial value of output records to match the value presently in the instrument
 - @init ‘exception handler’
 - Often the same as the read back protocol entry

```
getF {  
    out "\$1?";  
    in "%f";  
}  
setF {  
    @init { out "\$1?"; in "%f"; }  
    out "\$1 %f";  
}
```

Adding Stream Device Support

- Make changes to configure/RELEASE file
 - Add entries for streams and ASYN
 - IOCAPPS=/usr/local/iocapps/R3.14.12.3
 - ASYN=\$(IOCAPPS)/modules/soft/asyn/4-21-asd2
 - STREAMS=\$(IOCAPPS)/modules/soft/streamDevice/2-5-asd8

- Modify the application src/Makefile

.....

```
streams_DBD += base.dbd
```

```
streams_DBD += $(ASYN)/asyn.dbd
```

```
streams_DBD += $(ASYN)/drvAsynIPPort.dbd
```

```
streams_DBD += $(STREAMS)/dbd/stream.dbd
```

.....

```
streams_LIBS += asyn stream
```

.....

Adding Stream Device Support

- Make changes to application Db/Makefile
 - Add entries for the instruments and ASYN
 -
 - DB += StreamsExample.db
 - DB_INSTALLS += \$(TOP)/streamsApp/Db/streamEx.proto
 - DB_INSTALLS += \$(ASYN)/db/asynRecord.db
- This copies the database and the protocol file to TOP/db directory
- The protocol file must be stored in one of the directories listed in the environment variable STREAM_PROTOCOL_PATH

Adding Stream Device Support

- Modify the IOC startup script

```
epicsEnvSet ("STREAM_PROTOCOL_PATH", ":{TOP}/db")
```

```
.....
```

```
drvAsynIPPortConfigure("${USER}", "Device IP Address:Port", 0, 0, 0)
```

```
asynOctetSetInputEos("${USER}", -1, "Add Input Terminator Here")
```

```
asynOctetSetOutputEos("${USER}", -1, "Add Output Terminator Here")
```

```
## Load record instances
```

```
dbLoadRecords "db/myDatabase.db", "P=${USER}:,PORT=${USER},ADDR=0"
```

```
dbLoadRecords "db/asynRecord.db",
```

```
"P=${USER}:,R=device,PORT=${USER},ADDR=0,OMAX=10,IMAX=10"
```

```
.....
```

- P,R – PV name prefixes – PV names are \$(P)\$(R)name
- PORT – ASYN port name from corresponding devxxxConfigure command

Lab Session:

Control a Network Attached Device

- Host www.xxx.yyy.zzz – TCP Port 24742
- '\n' command terminator, '\r\n' reply terminator
- *IDN?
 - Returns device identification string (up to 100 characters)
- LOAD?
 - Returns three floating-point numbers separated by spaces (1, 5, 15 minute load average)
- VOLTS?
 - Returns most recent voltage setting
- CURR?
 - Returns current readback ($\pm 11A$)

Lab Session:

Control a Network Attached Device

- ON?
 - Returns the current on/off status
- ON [0,1]
 - Turns supply OFF/ON (0/1)
- VOLTS x.xxxx
 - Sets voltage ($\pm 10V$ range)