

Preparing the Cryoshipper for Shipping

Charging the Cryoshipper

Please note: when putting the lid on cryoshipper, make sure notch in the lid is lined up with sample holder handle.

Completely fill the cryoshipper with liquid nitrogen. Let it sit overnight. Top off with liquid nitrogen next morning. Repeat until no additional liquid nitrogen is absorbed.

Adding samples

Always use the shipper with the metal insert.

Shipping

Before shipping the cryoshipper, briefly remove the sample holder and **pour out all the excess liquid nitrogen**. You cannot ship the cryoshipper when it still contains liquid nitrogen in the central cavity. The phrase to use for couriers (e.g. FedEx) is “dry shipper”. The unit complies with all IATA regulations.