

# X-ray Diagnostics



Petr Ilinski  
Accelerator Division, NSLS-II, BNL  
NSLS-II ASAC  
February 1-2, 2012


# X-ray diagnostics objectives

---

- **e-beam**
  - emittance
    - source size measurements
      - pinhole
      - CRL
  - particles energy spread
- **Photon beam**
  - commissioning
    - ID radiation centering
  - position monitors
    - feedback

# e-beam Diagnostics

# e-beam Source Size



U.S. DEPARTMENT OF  
**ENERGY**

P. Ilinski, X-ray Diagnostics, NSLS-II ASAC, Feb. 1-2, 2012


4

**BROOKHAVEN**  
NATIONAL LABORATORY  
BROOKHAVEN SCIENCE ASSOCIATES


# e-beam / Radiation

	<b>BM</b>	<b>3-pole wiggler</b>
$\sigma_x [\mu\text{m}]$	75	140
$\sigma_y [\mu\text{m}]$	14	12
$\sigma_{x'} [\mu\text{rad}]$	56	77
$\sigma_{y'} [\mu\text{rad}]$	0.8	0.9
B [T]	0.4	1.14
$E_c [\text{keV}]$	2.4	6.8
$\sigma'_{ph} \sim 1/\gamma [\mu\text{rad}]$	170	170
$\sigma_{ph} [\mu\text{m}]$	0.24	0.09
$L_f [\text{mm}]$	4.3	1.5


# Radiation Sources - Spectral Flux



# Imaging with Pinhole



EFFECTIVE RESOLUTION IN MIXED MODEL



Pinhole diameter =  $2a$

$$\frac{1}{f} = \frac{1}{z} + \frac{1}{d} \quad \frac{1}{f} = \frac{(1+M)^2}{LM}$$

$$L = z + d \quad z = L/(1+M)$$

$$M = d/z \quad d = ML/(1+M) \quad M = \frac{L}{z} - 1$$

$$F = \frac{a^2}{\lambda f}$$

- For an optimal  $F=0.7$   
 $a = \sqrt{0.7 \lambda LM / (1+M)}$   
 and PSF (w=WFHM) at  
*Image Plane*:

$$w \approx 2 \sqrt{\frac{\lambda}{2\pi}} \frac{d^2}{f}$$

$$w \approx \frac{2}{\sqrt{2\pi}} \sqrt{\lambda M L}$$

*Source Plane*:

$$w \approx 0.7 \sqrt{\lambda L/M}$$

$$\sigma \approx 0.3 \sqrt{\lambda L/M}$$


$$\sigma \approx 0.3 \sqrt{\lambda Lz / (L-z)}$$




# Square Pinhole - PSF

E=12 keV, z=6m, d=10m

$5 \times 5 \mu\text{m}^2$




$10 \times 10 \mu\text{m}^2$




$33 \times 33 \mu\text{m}^2$

**Optimal**



$100 \times 100 \mu\text{m}^2$




U.S. DEPARTMENT OF  
**ENERGY**

P. Ilinski, X-ray Diagnostics, NSLS-II ASAC, Feb. 1-2, 2012 8


**BROOKHAVEN**  
NATIONAL LABORATORY  
BROOKHAVEN SCIENCE ASSOCIATES

# Pinhole setup optimization

## Resolution (PSF WFHM) vs. Pinhole position from the source



# X-ray Optics - Resolution



Larger NA allows to increase spatial resolution



A Compound Refractive Lens (CRL) optics will be used for high resolution setup

# Compound Refractive Lens (CRL) - PSF




## Calculated Point Spread Functions (PSF)

BM radiation, 0.4T, single electron,  $E_{ph}=20$  keV,  $\sigma'=0.05$  mrad

CRL - 12 Be lenses, 1-mm-diam,  $R=50\mu m$ ,  $z=3.02m$ ,  $d=12.98m$ ,  $M=4.3$



# X-ray Detection - PSF



## Spatial resolution is limited by:

- Diffraction:  $R \sim \lambda/NA$
- Depth of focus:  $R \sim \delta z \cdot NA$
- Spherical aberrations of substrate:  $R \sim t \cdot NA^3$
- Photoelectron range

Koch et al, J.Opt.Soc.Am. A, Vol. 15, p.1940, 1998



U.S. DEPARTMENT OF  
**ENERGY**

# X-ray Detection - beamline setup


---

- Required spatial resolution
  - 2-3  $\mu\text{m}$
- Field of View (FoV)
  - is defined by source size, required resolution, objective magnification, CCD dimensions
  - source is highly asymmetric
  - FoV needs to be 2÷3 times of the image size ( $6\sigma$ )
- FoV vs. Resolution
  - Different magnification for Horizontal and Vertical directions allows to achieve required resolution with standard CCD
  - Horizontal and Vertical slits
  - Linear CRL for imaging in Horizontal and Vertical planes


# Optical Layout - Resolution, FoV

	Units	BM	BM 3-DW	BM 8-DW	3PW	3PW 3-DW	3PW 8-DW
Source size V	um	14	14	14	12	12	12
Source size H	um	75	49	39	133	137	138
Aspect ratio H/V		5	4	3	11	11	12

Optical Element	Material	Delta	BL	Energy	D_eff	R	N	f	z	d	L= z +d	PSF	M= d/z	PSF /M	Source Size Hor	Source Size Vert	FoV= 6s*2*M Hor	FoV= 6s*2*M Vert	CCD FoV	CCD # pixels	CCD pixel size	Objective M	Objective Resolution	CCD pixel size/M
CRL				keV	um	um		m	m	m	um	um	rms	um	rms um	mm	mm	mm		um	x	um	um	
V-CRL	Be	8.5E-7	BM	20	420	50	12	2.45	3.02	12.98	16	1.0	4.3	0.2	75	14	3.9	0.7	0.9	1000	9.0	10	0.9	0.9
H-CRL	Be	8.5E-7	BM	20	428	50	8	3.68	10.28	5.72	16	0.6	0.6	1.1	75	14	0.5	0.1	0.9	1000	9.0	10	0.9	0.9
V-CRL	Be	8.5E-7	3PW	20	430	50	7	4.20	5.66	16.34	22	1.8	2.9	0.6	140	12	4.9	0.4	0.9	1000	9.0	10	0.9	0.9
H-CRL	Be	8.5E-7	3PW	20	428	50	8	3.68	17.33	4.67	22	0.5	0.3	1.9	140	12	0.5	0.0	0.9	1000	9.0	10	0.9	0.9
Slit																								
V-Slit			BM	20	20				2.85	13.15	16	20.1	4.6	4.5	75	14	4.2	0.8	1.8	1000	9.0	5	2.0	1.8
H-Slit			BM	20	26				6.00	10.00	16	12.4	1.7	7.3	75	14	1.5	0.3	1.8	1000	9.0	5	2.0	1.8
V-Slit			3PW	20	25				4.69	17.32	22	20.1	3.7	5.6	140	12	6.2	0.5	1.8	1000	9.0	5	2.0	1.8
H-Slit			3PW	20	31				10.00	12.00	22	12.0	1.2	10.0	140	12	2.0	0.2	1.8	1000	9.0	5	2.0	1.8




# Flux Budget




Source	I	H	V	Energy	BW	X-ray Flux	Window	Filter	X-ray Flux	Scintillator	Scint thickness	V-ph Flux	Objec tive	V-ph Flux	spot size H	spot size V	V-ph Flux/ pixel
	mA	urad	urad	keV	eV	ph/s			ph/s		um	ph/s	NA	ph/s	pixels	pixels	ph/s/px
BM 0.3T	500	10	10	all	all	1.4E+14	-	N2 10m	1.6E+12	CdWO4	100		0.3		300	300	
BM 0.3T	500	10	10	all	all	1.4E+14	Al 1mm	N2 20m	2.4E+11	CdWO4	100	2.6E+14	0.3	1.2E+12	300	300	1.3E+7
BM 0.3T	500	10	10	20	20		Al 1mm	N2 20m		CdWO4	100	4.4E+11	0.3	2.0E+9	300	300	2.3E+4
BM 0.3T	1	10	10	20	20		Al 1mm	N2 20m		CdWO4	100	8.8E+8	0.3	4.1E+6	300	300	4.5E+1
BM 0.3T	500	10	10	all	all		Al 1mm	N2 20m		CdWO4	10	9.3E+13	0.3	4.3E+11	300	300	4.8E+6
BM 0.3T	500	10	10	20	20		Al 1mm	N2 20m		CdWO4	10	1.2E+11	0.3	5.6E+8	300	300	6.2E+3
BM 0.3T	1	10	10	20	20		Al 1mm	N2 20m		CdWO4	10	2.4E+8	0.3	1.1E+6	300	300	1.2E+1
BM 0.4T	500	10	10	all	all	2.5E+14	-	N2 10m	1.0E+13	CdWO4	100		0.3		300	300	
BM 0.4T	500	10	10	all	all	2.5E+14	Al 1mm	-	5.0E+12	CdWO4	100		0.3		300	300	
BM 0.4T	500	10	10	all	all	2.5E+14	Al 1mm	N2 20m	2.5E+12	CdWO4	100	3.0E+15	0.3	1.4E+13	300	300	1.5E+8
BM 0.4T	1	10	10	all	all	5.0E+11	Al 1mm	N2 20m	5.0E+9	CdWO4	100	6.0E+12	0.3	2.8E+10	300	300	3.1E+5
BM 0.4T	500	10	10	20	20		Al 1mm	N2 20m	4.4E+9	CdWO4	100	5.5E+12	0.3	2.6E+10	300	300	2.8E+5
BM 0.4T	1	10	10	20	20		Al 1mm	N2 20m	8.8E+6	CdWO4	100	1.1E+10	0.3	5.1E+7	300	300	5.7E+2
3PW	500	10	10	all	all	2.0E+15	-	N2 20m	9.6E+14	CdWO4	100		0.3	4.6E+15	300	300	5.1E+10
3PW	500	10	10	20	20	1.3E+12	-	N2 20m	8.7E+11	CdWO4	100		0.3	3.3E+13	300	300	3.7E+8
3PW	1	10	10	20	200	2.6E+09	-	N2 20m	1.7E+09	CdWO4	100		0.3	6.6E+10	300	300	7.3E+5
3PW	500	10	10	all	all	2.0E+15	Al 1mm	N2 20m	6.0E+14	CdWO4	100	8.3E+16	0.3	3.9E+14	300	300	4.3E+9
3PW	500	10	10	20	20	1.3E+13	Al 1mm	N2 20m	3.9E+11	CdWO4	100	4.8E+14	0.3	2.2E+12	300	300	2.5E+7
3PW	1	10	10	20	20	2.6E+10	Al 1mm	N2 20m	7.8E+08	CdWO4	100	9.6E+11	0.3	4.5E+9	300	300	5.0E+4

# X-ray Diagnostics Beamlines Layout



# BM-A Layout (girder)




U.S. DEPARTMENT OF  
**ENERGY**

P. Ilinski, X-ray Diagnostics, NSLS-II ASAC, Feb. 1-2, 2012 17

**BROOKHAVEN**  
NATIONAL LABORATORY  
BROOKHAVEN SCIENCE ASSOCIATES

# BM-A Layout (girder)




U.S. DEPARTMENT OF  
**ENERGY**

P. Ilinski, X-ray Diagnostics, NSLS-II ASAC, Feb. 1-2, 2012 18

**BROOKHAVEN**  
NATIONAL LABORATORY  
BROOKHAVEN SCIENCE ASSOCIATES

# 3PW Layout (girder)




U.S. DEPARTMENT OF  
**ENERGY**

P. Ilinski, X-ray Diagnostics, NSLS-II ASAC, Feb. 1-2, 2012 19

**BROOKHAVEN**  
NATIONAL LABORATORY  
BROOKHAVEN SCIENCE ASSOCIATES

# 3PW Layout (girder)




U.S. DEPARTMENT OF  
**ENERGY**

P. Ilinski, X-ray Diagnostics, NSLS-II ASAC, Feb. 1-2, 2012 20


**BROOKHAVEN**  
NATIONAL LABORATORY  
BROOKHAVEN SCIENCE ASSOCIATES

# Insertion Devices Commissioning

# ID radiation - Dipole Chamber




Vertical cross-section



U.S. DEPARTMENT OF  
**ENERGY**

# Diamond Screens for NSLS-II commissioning




# Diamond Luminescent Screens

## Diamond X-ray beam monitors

Diamond Materials, Freiburg, Germany


Available options:

- Type Free-standing diamond fluorescence screen
- Options Defined sub-micron surface roughness to scatter the light from laser
- Material Optical transparent diamond, doped for increased fluorescence yield
- Thickness  $> 20 \mu\text{m}$
- Size Up to 75 mm; rectangular, circular or elliptical



# Diamond Luminescent Screens at PETRA-III

CVD Diamond Screens for PETRA III Beamlines  
M. Degenhardt, U. Hahn and H. Schulte-Schrepping  
HASYLAB Annual Report, DESY 2009




Alignment of the PETRA-III undulator beam in the pre-aligned front-end of beamline P09

- a) Bending magnet radiation, b) Bending magnet and undulator radiation,
- c) ongoing alignment, d) properly aligned undulator beam


# X-ray Beam Position Monitors

# Typical NSLS-II Front End Configuration



U.S. DEPARTMENT OF  
**ENERGY**


# Thermal Stability of XBPM Head



100 watts applied to pre-mask  
Common support structure  
(thermally coupled)

100 watts applied to pre-mask  
Separate support structures  
(thermally de-coupled)

# XBPM Assembly Layout




Invar slit stand	inches	mm	change degrees C	change degrees F	MATERIAL	in/in/deg F	EXPANSION inches	EXPANSION mm	EXPANSION microns
top and bottom plates	2.0	50.8	0.2	0.36	INVAR	7.2E-07	5.2E-07	1.3E-05	0.00
stand tube	35.3	895.4	0.2	0.36	INVAR	7.2E-07	9.2E-06	2.3E-04	0.23
stage	3.5	88.9	0.2	0.36	ALUMINUM 6061	1.3E-05	1.6E-05	4.2E-04	0.42
slit body	1.8	44.5	0.2	0.36	COPPER	9.8E-06	6.2E-06	1.6E-04	0.16
threaded rods	3.5	88.9	0.2	0.36	INVAR	7.2E-07	9.1E-07	2.3E-05	0.02
									0.83




# Power density and XBPM Signal distribution

NSLS-II IVU-20, 10m, K=1.8, Tungsten Blades

Incident Power density [W/mm<sup>2</sup>] (vs. x,z)




XBPM signal density distribution [ $\mu\text{A/mm}^2$ ] (vs. x,z)



# Blade Temperature Distribution

Blade Optimization  
Tungsten, 0.2mm FBM, 0.2mm BNL



# Conclusions

---

- **e-beam diagnostics beamlines**
  - research - completed
    - two interchangeable setups (pinhole/CRL)
  - beamline design - 80% completed
- **Diamond Luminescent Screens**
  - PETRA-III blueprints
- **X-ray BPMs**
  - optimization - ongoing
  - assembly - in design
  - FMB - in contact