



# Infrared Imaging and Microspectroscopy

---

**Lisa M. Miller**


*National Synchrotron Light Source*

*Brookhaven National Laboratory*

*Upton, NY 11973*





**U.S. DEPARTMENT OF ENERGY**  
**OFFICE OF BASIC ENERGY SCIENCES**




**Funded under contract: DE-AC02-98CH10886**

**BROOKHAVEN**  
NATIONAL LABORATORY  
BROOKHAVEN SCIENCE ASSOCIATES



# Chemical Features of Biological Components




## Nucleic Acid



## Carbohydrate




# Protein Structure Determination with FTIR



## Amide I Secondary Structure Assignments:

1620 - 1640	$\beta$ -sheet
1645	extended coil ( $D_2O$ )
1648 - 1657	$\alpha$ -helix
1660	triple helix
1670 - 1695	$\beta$ -sheet, $\beta$ -turn



**Collagen**  
(triple helix)

**Myoglobin**  
( $\alpha$ -helix)


**Abeta Amyloid**  
( $\alpha$ -helix +  $\beta$ -sheet)

**Cytochrome c**  
( $\alpha$ -helix + coil)


**STI**  
( $\beta$ -sheet + coil)

# IR Spectroscopy vs. Micro-Spectroscopy

## Bulk samples





## Microscopic Heterogeneity




# Why SYNCHROTRON Infrared Microspectroscopy?

Biological samples are **SMALL**:


- spatially
- in concentration



## Synchrotron Advantages:





# Single Red Blood Cell: Globar vs. Synchrotron Source



Diffraction Limit	
4000 cm <sup>-1</sup>	2.5 μm
2950 cm <sup>-1</sup>	3.4 μm
1650 cm <sup>-1</sup>	6.1 μm
1200 cm <sup>-1</sup>	8.3 μm
1000 cm <sup>-1</sup>	10 μm
600 cm <sup>-1</sup>	17 μm
200 cm <sup>-1</sup>	50 μm


- Synchrotron IRMS is diffraction-limited
- Conventional IRMS is throughput-limited








# Beamline U10B

- Synchrotron infrared beam pipe
- Nicolet Magna 860 FTIR
- Spectra Tech Continuum IR microscope
- Filter cube turret
- UV (quartz) light source
- Schwarzchild 32x IR objective
- Automated mapping stage



# Chemical Content of New Bone in Osteoporosis




- Rate of bone mineral formation is slower in osteoporosis






L.M. Miller, J. Tibrewala, C.S. Carlson. *Cell. Mol. Biol.*, **46**:1035-44 (2000).



# Prion Protein Structure and Location



- Hamster scrapie 263K
- Dorsal root ganglia
- Terminally infected



# Chemical Imaging of Hair




cuticle (3-5  $\mu\text{m}$ )




cortex (40-100  $\mu\text{m}$ )

medulla (5-10  $\mu\text{m}$ )



## Phospholipids





## Proteins



# Protein Aggregate Formation in Alzheimer's Plaques





- misfolded, aggregated protein is associated with Abeta plaques





L.M. Miller, P. Dumas, N. Jamin, J.-L. Teillaud, J. Miklossy, L. Forro (2002). *Rev. Sci. Instr.*, **73**: 1357-60.









# Metal Accumulation in Alzheimer's Plaques



- Zn and Cu are associated with Abeta plaques





# Correlation Between Metal Accumulation and Protein Aggregation






# Chemical Imaging of Mitotic Cells


OPTICAL IMAGE



PROTEIN ( $1650\text{ cm}^{-1}$ )








LIPIDS ( $2925\text{ cm}^{-1}$ )



N. Jamin, J.L. Teillaud, P. Dumas, G.L. Carr, G.P. Williams (1998). *PNAS* **95**: 4837-40.



# Chemical Changes in Apoptotic Cells




How does apoptosis proceed in real time?

L.M. Miller, P. Dumas, N. Jamin, J.-L. Teillaud, J. Miklossy, L. Forro (2002). *Rev. Sci. Instr.*, **73**: 1357-60.

# Impact of NSLS-II: Schematic Brightness Comparison

IR flux is primarily dependent upon ring current in mid-IR



- Mid IR: NSLS-II will have equal brightness to NSLS
- Far-IR and multi-point imaging: disadvantaged by NSLS-II

# Summary

## Overview

- Infrared microspectroscopy and imaging can be used to image the chemical makeup of biological tissues and cells
- The spatial resolution of synchrotron IRMS is diffraction-limited at 2-20  $\mu\text{m}$

## Importance of Research

- Chemical imaging of mineralized tissues, *in situ* protein, lipid, nucleic acid content, protein structure

## Scientific Challenges and Opportunities

- Multi-technique imaging: combining x-ray, visible, IR imaging techniques
- Improve spatial resolution
- Improve data collection rates

## Impact of NSLS-II

- NSLS-II brightness will slightly improve the mid-IR imaging capabilities
- NSLS-II will prevent the creation of large opening-angle beamlines for (1) far-IR applications, and (2) line-source (i.e. multi-point) imaging

## Synergy and Demand

- Combination of IRMS with x-ray imaging techniques
- Currently 4 IRMS bending magnet beamlines; need at least this many for NSLS-II