

Scientific Opportunities with Photoemission Spectroscopy: Present and Future

Zahid Hussain

Division Deputy for Scientific Support

Advanced Light Source

Lawrence Berkeley National Laboratory



Outline



- The big picture
- The leading edge of measurements today
 - some examples
- The future:
 - Spin-Resolved photoemission using TOF analyzer and exchange scattering
 - -nm-scale measurements, "nanoARPES"
 - -time-resolved photoemission
- · The MERLIN beamline

Understanding complex correlated phenomena require sharper and sharper tools



The Big Picture

Grand Challenges & Why x-rays

Grand Challenges!!



Energy problem - search for 20 TWatts of energy,

solar energy, hydrogen fuel, nuclear energy ??

Membrane Proteins – from 3D structure of
Macromolecules to understanding functions-dynamics
(engineering of drugs)

Understanding Emergent Phenomena -

Nanoscale materials, Strongly correlated electron systems (high Tc superconductor.....)


The ultra-small - imaging of single atom in a host of others, manipulation of single spin + ...

The ultra-fast - science in the ps-fs-as

Why X-Rays (& not neutrons or electrons)




Tunable x-rays offer variable interaction cross section



Science with Light Sources





Adopted from: Franz




Emergent Phenomena

Strongly Correlated electron Systems

Strongly Correlated Electronic Systems

Condensed matter physics







Colossal Magnetoresistance (CMR) Effect



Novel Electronic Phases





CO: Charge Order (Stripes)

FI: Ferromagnetic Insulator

AFI: Antiferro. Insulator

CAF: Canted AFM Insulator


CMR: Colossal MagnetoResis.

- O Large drop of resistivity upon relatively small magnetic fields
- O Para → Ferromagnetism
- O Most dramatic on the insulating phase (short range orbital order)

Manganites Exhibit Interplay of Charge, Spin, lattice and Orbital degrees of freedom



Interacting degrees of freedom (complex electron systems)




Competition among many <u>Energy</u> and <u>Length</u> scales

Determine the physics of these systems

Energy Scale of Important Excitations





- >Superconducting gap ~ 1 100meV
- >Optical Phonons: ~ 40 70 meV
- > Magnons: ~ 10 meV 40 meV
- > Pseudogap ~ 30-300 meV
- >Multiphonons and multimagnons ~ 50-500 meV
- > Orbital fluctuations (originated from optically forbidden d-d excitations): ~ 100 meV 1.5 eV

Requirement: High Energy Resolution and High Efficiency Detection

Fundamental Spectroscopies of Condensed Matter



Spectral functions (One-particle properties)
Correlation functions (two-particle properties)

1-particle response


•Angle resolved photoemission (ARPES) : Single-particle spectrum $A(k,\omega)$




• <u>Spin</u>: Inelastic Neutron Scattering (INS): (neutrons carry magnetic moment)
Spin fluctuation spectrum $S(q, \omega)$

• Charge : Inelastic x-ray scattering (IXS) : Coupled excitation in the Charge Channel $N(q, \omega)$

(MERLIN/QERLIN (ALS); FEL)







?

Angle-resolved Photoemission Spectroscopy










Scientific issues

What does ARPES measure?







Properties of the photoelectron reflect the properties of the injected hole

namely,


we can learn about the scattering and lifetime of the injected hole

Courtesy: Eli Rotenberg

"kinkology"



The carriers have a finite lifetime due to absorption and emission of phonons and other excitations




[1] Koralek et al., Phys. Rev. Lett. 96, 017005 (2006)

Quasiparticle in Manganite: LSMO x=0.4




- Nodal quasiparticle:
- nodal-antinodal dichotomy in a *non* superconducting material



· Are the pseudogap state & the -0.4 0 nodal-antinodal dichotomy hallmark of the

superconductivity state?




N. Mannella et al. Nature <u>438</u>, 474 (2005)

LAWRENCE BERKELEY NATIONAL

Temperature evolution of the small QP peak linked to transport properties




and the metal-insulator transition in LSMO Mannella et al



How could measurement of a microscopic electronic structure in certain part of the BZ could be related to the bulk macroscopic property?

Spin-Charge Separation





Photohole decay into

Two topological defects




Spinon: speed controlled controlled

holon: speed


by t (nearest-neighbor hopping

(exchange Coupling constant) Amplitude)



Spin-Charge Separation in 1D SrCuO₂






Holons and spinons have different mass & thus different dispersion (Fermi velocity).

B.J. Kim, E. Rotenberg' C. Kim et al. Nature Physics 2, 355(2006)

Band Structure of C_{60} on Ag (111) and Ag (100) Surfaces - Orientation Dependence






Same Hexagon Structure; Completely Different Dispersion!

Brouet, Yang, Hussain, Shen et al, Science, <u>300</u>, 303 (2003); PRL (2004)

Combination of Experiment and Theory (strong orientation dependence)







Dotted Lines: Theory (Louie, Cohen et al) LDA);

2D Images: Experiments

C60 (Conclusion)








Wang, Brouet, Louie, Hussain, Shen et al Science 300, 303(2003) PRL93, (2004)


AWRENCE BERKELEY NATIONAL LABORATORY


Purely Two Dimensional State: Graphene



Graphene: A single hexagonal sheet of Carbon atoms that is the building block of carbon nanotubes, graphite, buckyballs, and other C-based materials







A gapless semiconductor with a famous conical bandstructure.

- ballistic transport at room temperature
- high current capacity (130 nA per nanotube, or 100,000,000 A/cm²


By imposing a gate voltage on a graphene layer we can easily tune the doping from *n* to *p* type, and we can realize a new generation of ultra-small, high current capacity devices.

LABORATORY

Example 2: Breaking down the barriers to carbon-based electronics




- The electronic properties of graphene depend on the thickness and applied field.
- Demonstrating the possibility of 0.6 nm thick switches based on carbon.




Ohta, T., Bostwick, B., Seyller, T., Horn, K. & Rotenberg, E. Controlling the Electronic Structure of Bilayer Graphene. Science 313, 951-954 (2006).

Photoexcited Electron Emission from Diamond and Diamondoid-thiol Self-assembled Monolayer





Himpsel etc. 1979 (Wisconsin) Pate etc. 1986 (SSRL)



Tetramantanethiol (4-cage)


- Self Assemble on Metal Surface (Ag, Au, ...)
- Natural resource (oil field)
- More uniform
- 5 times more effective
- Better electron conductivity

WL Yang, Z Hussain, ZX Shen, Science 2007 (ALS)

Patent filed, Stanford-LBNL-Chevron


Surface Destruction Test





Electron Emission from Diamondoid






Space Charge Effect in Photoemission

Caution





X. J. Zhou, B. Wannberg, W. L. Yang, V. Brouet, Z. Sun, J. F. Douglas, D. Dessau, Z. Hussain, Z.-X. Shen, J. Electron Spectroscopy and Related Phenomena, 142, 25 (2005).




What are new opportunities with core level photoemission?


In-Situ & dynamical studies of chemical reactions at surfaces

BERKELEY LAB

Photoelectron Spectroscopy








ag.zh/photospectros1/7-97

High resolution C1s Photoelectron Spectra of hydrocarbons (sharper tool with modern SR)





C 1s photoelectron spectra of propyne (res: SR~30meV, lab source ~ 0.5 eV (dashed))

Unambiguous assignment of peaks in propyne spectrum is made possible by characteristic vibrational structure of model compounds (ethyne and ethane) and abinitio theory.

Shift of the methyl (CH3) peak in propyne relative to ethane is due to the electronegativity of the ethyne (HC°C) group.


From BL 10.0 (AMO, ALS) Thomas et al, PRL

Ambient pressure photoemission

Ionization Energy (eV)
LAWRENCE BERKELEY NATIONAL LABORATORY

Ambient Pressure Photoemission: ~ 10 torr






Oxidation dynamics



Strong temperature dependence




Si(100) oxidized by water vapour @ .1 torr

(S. Mun et al)

Steps Toward Hydrogen Vehicles





Nikitin et. al., *Phys Rev Lett.* **95**, 225507 (2005)

Characterising new materials:

Carbon Nanotubes for storing hydrogen: safely, efficiently and compactly.

The DOE Freedom CAR program has set the goal of a material that can hold 6% of the total weight in hydrogen by the year 2010. Theoretical calculations indicate they may exceed these goals substantially.





Funded by DOE, NSF and Global Climate and Energy Project (alliance of scientific researchers and leading companies in the private sector, including ExxonMobil, General Electric and Schlumberger)

The Ultimate Solution: Energy Problem (Hope!)



Need study of transition metal oxides as Photocatalysts for breaking the water using sun light (tailor properties using emergent phenomena and quantum confinement.








■TiO₂ as a Energy Source

Hydrogen Production from the Water ($(H_2O) + photons \rightarrow H_2 + \frac{1}{2}O_2$)

(problem: only absorbs UV light, look for other catalyst)


- High payback!
 - •Unlimited source (water + solar light)
 - Very Clean and Recyclable Exhausts (Fuel Cell)



courtesy:Akira

Fuel Cell - Schematic






Anode: Hydrogen oxidation

Hydrogen gas = Hydrogen Ions + Electrons






Cathode: Oxygen reduction
Oxygen gas + Protons + Electrons = Water

In both cathode and anode, Pt based catalysts are applied to increase the rate of each chemical reactions. Need better material than presently used Pt.

Cathode: the performance of polymer electrolyte membrane fuel cells is limited by the slow rate of O_2 reduction (ORR) at Cathode, ~5 orders of magnitude slower than H_2 oxidation at Anode

Breakthrough Research For Fuel Cells





Hydrogen powered fuel cells for automotives are one of the cleanest and most efficient technologies for generating electricity. But until now, progress is slow.

The discovery of a unique <u>platinum-nickel</u> <u>alloy</u> represents a breakthrough in catalyst research: it is <u>90 times more active</u> than state-of-the-art platinum catalysts currently used.

Science 315 Jan 2007

In standard Pt catalysts (top) absorption of oxygen on the surface is hindered by the binding of other molecules, such as OH. In the new material (bottom) The nickel atoms change the surface properties such that OH cannot bind as well, leaving room for oxygen.

Research team includes: Argonne and Berkeley National Labs, U. South Carolina.







Scientific opportunities with ARPES

for unraveling of complex problems

X-ray Spectroscopy of Condensed Matter



Quantum Number Selectivity:

✓ Absorption

$$\omega \, \varepsilon_2 \Rightarrow \Delta E = E_f - E_{\iota}$$

✓ Angle-integrated photoemission

$$N(E,h\omega) \Rightarrow E_{f,E_{\iota}}$$

✓ Angle-resolved photoemission (also inelastic scattering)

$$N(E, h\omega, \theta, \varphi) \Rightarrow E_{f, E_{\iota;}} \kappa$$

!!! Spin-polarized photoemission


$$(N_{\uparrow} - N_{\downarrow})/(N_{\uparrow} + N_{\downarrow}) \Rightarrow E_{f,E_{\iota,}} \kappa, \sigma$$

Photoemission with circularly polarized light and spin detection



Selective excitations (use of elliptically polarizing undulator)





Courtesy: Yulin Chen

Spin detection (two schemes)



Mott Detector Spin-orbit interaction

Hint = $L \cdot S$

e- counter

e- 20-100 KeV

Gold target
(heavy nuclei)

FOM ≤ 10⁻⁴


Exchange scattering interaction

Reflectivity contains a

term:

 $\alpha P \cdot M$

e- counter



Magnetized thin film

(e.g.,Co/W(110))

x 100

FOM ~ 10⁻²

R. Bertacco et al. 2001 Hillebrecht et al. 2002

R. Zdyb and E. Bauer 2003

D.T. Pierce et al. 1988 +....

Spin-Resolved Photoemission (TOF Project)



"Time-of-Flight" energy analysis

Multichannel detection in time (energy):

~ 10-100 times more efficient than single channel dispersive analyzer

"Exchange Scattering" based spin analysis

- ~ 100 times more efficient than Mott Detector
- > Spin-Resolved ARPES: Improved efficiency & high resolution
 - > Overall FOM: ~ 1000 times vs. existing (Mott det.+ dispersive analyzer)

energy resolution : ~ 10meV


Better statistics : > factor of 10

Note: TOF is inherently low noise detection as detector counts for short time; only the time window when electrons of interest arrive craf, Schmid, Jozwiak, Hussain,

Lanzara LAWRENCE BERKELEY NATIONAL LABORATORY et al, PRB, 71, 144429 (2005)

Spin asymmetry (A_s) : Co/W(110) Experimental results using SPLEEM







Graf, Schmid, Jozwiak, Lanzara, Hussain et al, PRB, 71, 144429 (2005)

TOF Spin-Resolved Photoemission








Nano ARPES






(lectured by Janos Kirz)


Time-Resolved Photoemission Comparison of the Hemispherical Analyzer and the TOF Analyzer proposed







Currently used Hemispherical Analyzer (2D detection)

(Bi2212 Bi-layer splitting)







Key Specifications:
Acceptance angle: 30 degree

Energy resolution: <=2meV (5eV Pass Energy)

Angular resolution: <=0.1 degree (~2mrad)

(comparable to Scienta analyzer but 100 times faster)


Fabrication – in progress

Time-Scale of Various Phenomena



- Ultra-fast time regime:
 < 200fs (ΔE >10meV)
 - Electron excitation/de-excitation (fs)
 - Bond breaking
 - · Carrier-carrier scattering
 - Hole-optical phonon scattering
 - Charge density wave/charge transfer
 - Magnetic Dynamics
 - Relaxation of biological system after light absorption (Rhodopsin);
- Time regime: ~ 200fs 2ps (∆E ~1-10 meV)
 - Phase transition (diamond → graphite)
 - Carrier acoustic phonon scattering
- Time regime: > 1-100ps
 - Stripe fluctuation in High Temp Superconductor
 - Magnetic recording
 - Protein folding (ps-s)

meV Resolution Soft X-Ray Beamline (MERLIN)




Specifications:

- >Resolving power: E/ΔE=100,000 with 5μm slits
 - i.e. better than 1 meV when photon energy is below 100eV
- >Photon energy range: 15eV to 100eV, fully optimized
 - maximum achievable photon energy ~140eV
- >Elliptically Polarized Undulador (EPU): full polarization selection (linear and/or circular)
- **Photon Flux:** ~ 5 ×10¹¹ photons/s/meV

Optical Layout (SGM)

- ✓Ultra High Resolution& High Resolution Modes
- ✓ Asymmetric Undulator to suppress higher orders



meV Resolution Spectroscopy (MERLIN) Beamline




Specifications:

- >Resolving power: E/ΔE=100,000 with 5μm slits
 - i.e. better than 1 meV when photon energy is below 100eV
- **≻Photon energy range:** ∼ 8eV to 100eV, fully optimized
 - maximum achievable photon energy ~150eV
- >Elliptically Polarized Undulador (EPU): full polarization selection (linear and/or circular)
- **Photon Flux:** ~ 5 ×10¹¹ photons/s/meV


MERLIN Layout

Ultra-High Resolution & High Resolution Modes
Quasi-periodic undulator for Higher order




Reduced Flux in Higher Harmonics





Ray Tracing (Shadow)





hv ~ 30eV

Groove density=1,800 l/mm

 ΔE_{slit} = 0.029meV/ μ m ΔE_{slope} = 0.22meV/ μ rad


$\Delta E_{\text{slit}} \text{=} 0.14 \text{meV} @ 5~\mu\text{m}$ slit

slope error (μrad)	0.25	0.50
ΔE _{slope} (meV)	0.055	0.11
$\frac{\Delta E_{\text{tot}}}{\Delta E_{\text{slope}}^2)^{1/2}} = (\Delta E_{\text{slit}}^2)^{1/2}$	+ 0.15	0.18

Yi-De Chuang

Resolution Rowland Mode






LEG: 900 l/mm

Defocus, coma = zero. Line curvature negligible.

LEG: 900 l/mm, 10µm slits MEG: 1800 l/mm 5µm slits HEG: 3600 l/mm, 5µm slit

Flux @ sample






ID flux into 0.9×0.9 mrad²
Transmission entrance slit
Reflectivity all mirrors
Grating efficiencies
Band width correction

Cutting-edge technology spherical grating substrate




High resolution grating substrate (measured at BESSY metrology laboratory)

 $0.23 \mu rad RMS slope error$



MERLIN R=15m, center line, profile of residual slope




residual slope error: 1.03 µrad pv / 0.23 µrad rms,

radius of curvature: 14.965 m

meV Resolution VLS Spectrograph








Ray Traces

- Calculated/measured Resolution
 3 meV (high efficiency)
- · Overall length = 2 meters.
- Spectrograph for Merlin beamline (completion summer 2007)





 $hv = 49 \text{ eV} \pm 5 \text{meV}$

Spectrograph: Energy Resolution Test



Straight beam with ~ 6µm source size BL 12.0






Low-Temperature Goniometer with Six Degrees of Freedom



Sample Position



Cryostat

Support tube

- (1). Six degrees of freedom: 3 rotational and 3 translational;
- (2). Samp temperature ~10 K (no radiation shield);
- (3). Stability of sample against temperature change

Designed and fabricated by John Pepper (BL 10.0)

AWRENCE BERKELEY NATIONAL LABORATORY

Acknowledgement



- Photoemission: Z X Shen (Stanford), Eli Rotenberg, (LBNL, ALS), Xingjiang Zhou (ALS, China), Wanli Yang (ALS), Norman Mannella (ALS, Stanford), Simon Mun (ALS), Akira (NSLS)
- o TOF Analyzer: Chris Jozwiac, Ale Lanzara, (MSD, UCB), Gennadi Lebedev, ZX Shen, Yulin Chen
- o Inelastic Scattering: Z. Hasan (Princeton), Yi-De Chuang, Jinghua Guo; Jonathan Denlinger, (LBNL, ALS), Eric Guillikson,
- o Nanoscience characterisation: Franz Himpsel (univ of Wisconsin),

Lou Terminello (LLNL), Andreas Scholl (ALS), Jo Stohr (SSRL)

o MERLIN Beamline: Malcolm Howells, Rubin Reninger, Nicholas Kellog, Yi-De Chuang, Z. Hassan, Allesandra Lanzara, ZX Shen