

Basics of Data Processing


Scott Calvin
Sarah Lawrence College

September 30, 2005

2005 NSLS EXAFS Course:
Theory, Experiment, and Advanced Applications

1

Importing Raw Data in Athena



September 30, 2005

2005 NSLS EXAFS Course:
Theory, Experiment, and Advanced Applications

2

What This Talk Is

- An overview of some of the general features of Artemis and Athena
- A highlight of a few areas where novices often go astray

September 30, 2005

2005 NSLS EXAFS Course:
Theory, Experiment, and Advanced Applications

3

What This Talk Is Not


- An exhaustive documentation of every feature in Artemis and Athena
- A detailed tutorial
- A discussion of EXAFS theory

September 30, 2005

2005 NSLS EXAFS Course:
Theory, Experiment, and Advanced Applications

4

Groups in Athena




"Groups" include individual scans

September 30, 2005

2005 NSLS EXAFS Course:
Theory, Experiment, and Advanced Applications

5

"Marked" and "Current" Groups



This is the "current" group

Note that these buttons plot the current group


These buttons plot all the marked groups

September 30, 2005

2005 NSLS EXAFS Course:
Theory, Experiment, and Advanced Applications

6

"Marked" and "Current" Groups




September 30, 2005

2005 NSLS EXAFS Course:
Theory, Experiment, and Advanced Applications

7

Compare Scans



Vertical translation between
scans will not affect EXAFS


But a closer look...

September 30, 2005

2005 NSLS EXAFS Course:
Theory, Experiment, and Advanced Applications

8

Compare Scans



...reveals other problems.

Energy calibration appears to have drifted


Scans 38 and 39 appear inconsistent with the rest

September 30, 2005

2005 NSLS EXAFS Course:
Theory, Experiment, and Advanced Applications

9

Alignment in Athena




Under the "align" menu, you can have Athena attempt to shift scans for you, or you can do it "by hand"

September 30, 2005

2005 NSLS EXAFS Course:
Theory, Experiment, and Advanced Applications

10


Alignment Matters



The scans are now aligned

If we had not aligned them, the average would show considerable differences:

Note missing shoulder on unaligned version



September 30, 2005

2005 NSLS EXAFS Course:
Theory, Experiment, and Advanced Applications

11

alignment is **very important**

September 30, 2005


2005 NSLS EXAFS Course:
Theory, Experiment, and Advanced Applications

12

Background Removal in Athena

This information applies to the current scan

Right clicking brings up a menu that lets you apply these settings to all scans, or to all *marked* scans




September 30, 2005


2005 NSLS EXAFS Course:
Theory, Experiment, and Advanced Applications

13


Choosing E_0 in Athena



Click here to see a plot of the derivative spectrum



Choose E_0 at first peak in the derivative spectrum



September 30, 2005

2005 NSLS EXAFS Course:
Theory, Experiment, and Advanced Applications

14

Choosing E_0 carefully is of **low importance** for EXAFS, as long as it is done consistently for all scans

Normalizing in Athena

Background removal

E0: 9651.998 ☒ E0 shift: 0 Rbkg: 1.0 ☒

Standard: None Background: Autobk

k-weight: 1 ☐ Edge step: 1.3449 ☒ fix step: ☒ flatten ☐

Pre-edge range: -150 ☒ to -75 ☒

Normalization range: 150 ☒ to 1252.349 ☒


Spline range: k: 0.5 ☒ to 10.040 ☒

E: 0.952 ☒ to 1352.337 ☒


Spline clamps: low: None high: Strong

Choose pre-edge and
normalization ranges

Good and Bad Normalization



These **pre-** and **post-** edge lines look OK




This **post-** edge line goes too high near edge

September 30, 2005


2005 NSLS EXAFS Course:
Theory, Experiment, and Advanced Applications

17

Merging Data in Athena



Merge aligned and normalized scans here



September 30, 2005

2005 NSLS EXAFS Course:
Theory, Experiment, and Advanced Applications

18

Normalization may be *very important* in some cases, such as determining coordination number when you think you know S_o^2 , but if you are fitting S_o^2 it is *relatively unimportant*

R_{bkg}

Background removal

E0: 9651.998 E0 shift: 0 Rbkg: 1.0

Standard: None Background: Auto

k-weight: 1 Edge step: 1.3449 fix step: flatten

Pre-edge range: -150 to -75


Normalization range: 150 to 1252.349

Spline range: k: 0.5 to 18.840

E: 0.952 to 1352.337

Spline clamps: low: None high: Strong

R_{bkg} controls how far up in the Fourier transform the background will try to suppress; i.e. how "wiggly" the background can be

R_{bkg}


September 30, 2005

2005 NSLS EXAFS Course:
Theory, Experiment, and Advanced Applications

21

 R_{bkg} is fairly important


Ideally, small changes in R_{bkg} should not significantly change the parameters you find via your fits

September 30, 2005

2005 NSLS EXAFS Course:
Theory, Experiment, and Advanced Applications

22

Fourier Transforms in Athena




September 30, 2005

2005 NSLS EXAFS Course:
Theory, Experiment, and Advanced Applications

23

Choosing k_{\max}



This data appears consistent for all scans. If there is trouble at high- k , it is systematic and reproducible.

September 30, 2005

2005 NSLS EXAFS Course:
Theory, Experiment, and Advanced Applications

24

Choice of k_{\max} should be of **low importance**

Ideally, small changes in choice of k_{\max}
should not change parameters you find
via your fits


September 30, 2005

2005 NSLS EXAFS Course:
Theory, Experiment, and Advanced Applications

25

Choosing k_{\min}

Choose k_{\min} where background
stops being strongly dependent on
small changes in background
parameters



September 30, 2005

2005 NSLS EXAFS Course:
Theory, Experiment, and Advanced Applications

26

Choice of k_{\min} should be of **low importance**


Ideally, small changes in choice of k_{\min} should not change parameters you find via your fits. In practice, fit is often **more sensitive** to k_{\min} than k_{\max}

September 30, 2005


2005 NSLS EXAFS Course:
Theory, Experiment, and Advanced Applications

27


k-weight




Unweighted



k-weight 2



k-weight 1



k-weight 3

September 30, 2005

2005 NSLS EXAFS Course:
Theory, Experiment, and Advanced Applications

28

Athena Reminder!

This affects the Fourier transform

This affects the k-space plot

Forward Fourier transform

k-weight: 1 dk: 0 window type: kaiser-bessel

Plotting parameters

plot multiplier: 1 y-axis offset: 0

Displaying parameters for group "merged data" ... done.

Data groups

- merge unshifted
- merge shifted
- r_bkg = 0.2
- r_bkg = 1.0
- r_bkg = 2.0
- r_bkg = 0.2 bkg
- r_bkg = 1.0 bkg
- r_bkg = 2.0 bkg
- merged data

Plot current group in

Plot marked group in

Plotting options


- chi*k*kw
- chi
- chi*k
- chi*k*2
- chi*k*3

September 30, 2005

2005 NSLS EXAFS Course:
Theory, Experiment, and Advanced Applications

29

Effect of k-weight on Fourier Transform



September 30, 2005

2005 NSLS EXAFS Course:
Theory, Experiment, and Advanced Applications

30

k-weight *should* be of **low importance**
 Ideally, different k-weights should not
 change parameters you find via your fits.

September 30, 2005

2005 NSLS EXAFS Course:
Theory, Experiment, and Advanced Applications

31

Entering a Crystal Model

Lattice parameters are given in Angstroms

This is the maximum distance from the absorbing atom (in Angstroms) that you want feff to calculate. Make it a few Angstroms larger than the largest path you plan to use, so that feff can calculate potentials properly

Get here from the Theory menu

Artemis gives coordinates of atoms in fractional coordinates; i.e. multiply by lattice parameters to get Angstroms

Atoms	feff.inp	Interpretation	Core	El	X	Y	Z	Ts
ZnO								
A	3.25010							
B	3.25010							
C	5.20710							
Alpha	90.00000							
Beta	90.00000							
Gamma	120.00000							
Unit size	6.00000							
Edge	K							
Shift vector	0.00000							
	0.00000							
	0.00000							

September 30, 2005

2005 NSLS EXAFS Course:
Theory, Experiment, and Advanced Applications

32

What if the Sample Isn't Crystalline?

- Short-Range Amorphous (e.g. glasses): Use a crystal structure which is locally similar on the Atoms page, but only use the inner paths for analysis
- Macromolecular: Use the space group P 1 and lattice parameters $A=B=C=100$. Then enter coordinates of atoms in units of hundredths of an Angstrom.

September 30, 2005

2005 NSLS EXAFS Course:
Theory, Experiment, and Advanced Applications

33

Atoms Generates a feff.inp File

```

TITLE ZnO
HOLE 1 1.0 * Zn K edge (9659.0 eV), second number is
* mphase,mpath,mfeff,mchi
CONTROL 1 1 1 1
PRINT 1 0 0 0
RMAX 6.0
*CRITERIA curved plane
*DEBYE temp debye-temp
*NLEG 8
POTENTIALS
* ipot z element
0 30 Zn
1 30 Zn
2 8 O
ATOMS
* x y z ipot tag dis
0.00000 0.00000 0.00000 0 Zn1 0.0
1.87643 0.00003 -0.56393 2 O1 1.6
-0.93824 1.62508 -0.56393 2 O1_2 1.6
0.00000 0.00000 -3.18745 2 O1_3 3.1
1.87643 0.00003 2.60355 1 Zn1_1 3.2
-0.93824 -1.62502 2.60355 1 Zn1_1 3.2
  
```

This file is generated automatically when you run Atoms, but you should look at it to see if it makes sense

This list of potentials can be changed, but the absorber must be 0 and you cannot skip numbers

These coordinates and distances are given in Angstroms; you should always check that first-shell distances and coordination numbers make sense to you

September 30, 2005


2005 NSLS EXAFS Course:
Theory, Experiment, and Advanced Applications

34

The Artemis Main Screen

The left part of the screen changes depending on what is highlighted in... the "Data & Paths" panel.

The right side of the screen mostly controls the plots... and the big green Fit button makes it go!



September 30, 2005

2005 NSLS EXAFS Course:
Theory, Experiment, and Advanced Applications

35

Feff Screen

Interpretation shows you a list of paths:

Black non-italicised paths are included in the fit


Paths with the text in brown are not currently included in the fit

Italicised black paths are not currently in the path list at all

Paths with brown backgrounds are direct scattering

Paths with blue backgrounds are focused

Right-clicking allows you to add paths to the list, include them in the fit, and more...




September 30, 2005

2005 NSLS EXAFS Course:
Theory, Experiment, and Advanced Applications

36

Data Screen

Clicking on a data set brings up a screen where you can choose k- and R-ranges for fit, k-weights, etc..



September 30, 2005

2005 NSLS EXAFS Course:
Theory, Experiment, and Advanced Applications

37

Guess, Def, Set Screen

This screen allows you to define various kinds of parameters:

Guess parameters are optimized during fitting

Def parameters are calculated during fitting


Set parameters are calculated before fitting

After parameters are calculated after fitting

Skip parameters are not used

Restraints force a parameter to be close to a specified value

"reff" is a special parameter equal to half of the nominal path length




September 30, 2005

2005 NSLS EXAFS Course:
Theory, Experiment, and Advanced Applications

38

Guess, Def, Set Screen

Clicking on a path lets you specify path parameters in terms of Guess, Def, Set parameters




September 30, 2005

2005 NSLS EXAFS Course:
Theory, Experiment, and Advanced Applications

39

Log Files Screen

Clicking on Fit in the Data & Paths screen brings up a screen that lets you view results from previous fits and compare them graphically



September 30, 2005

2005 NSLS EXAFS Course:
Theory, Experiment, and Advanced Applications

40