

# X-RAY FLUORESCENCE MICROPROBE (XFM)

Project Team: P. Bertsch<sup>1</sup>, E. Dooryhee<sup>2</sup>, G. Flynn<sup>3</sup>, J. Flynn<sup>4</sup>, M. Ginder-Vogel<sup>5</sup>,

A. Lanzirotti<sup>6</sup>, M. Marcus<sup>7</sup>, L. Miller<sup>2</sup>, P. Northrup<sup>8</sup>, T. Punshon<sup>9</sup>, M. Rivers<sup>6</sup>,


B. D.P. Siddons<sup>2</sup>, G. Smith<sup>10</sup>, S. Sutton<sup>6</sup>, R. Tappero<sup>2</sup>, T. Tokunaga<sup>7</sup>

<sup>1</sup>Univ. Kentucky, <sup>2</sup>Brookhaven National Lab, <sup>3</sup>SUNY Plattsburgh, <sup>4</sup>American Museum of Natural History, <sup>5</sup>Calera Corp., <sup>6</sup>Univ. Chicago., <sup>7</sup>Lawrence Berkeley National Lab, <sup>8</sup>Stony Brook Univ., <sup>9</sup>Dartmouth College, <sup>10</sup>Indianapolis Museum of Art







## TECHNIQUES AND CAPABILITIES

- XFM is an optimized three-pole wiggler beamline for the characterization of materials in an "as-is" state that are chemically heterogeneous at the micrometer scale via synchrotron induced X-ray fluorescence.
- XFM includes instrumentation for microbeam X-ray fluorescence ( $\mu$ XRF), diffraction ( $\mu$ XRD) and fluorescence computed microtomography (FCMT). However, it is optimized to provide users state-of-the-art microfocused Extended X-ray Absorption Fine Structure ( $\mu$ EXAFS) spectroscopy between 4 to 20 keV.
- XFM will trade-off beam size and flux for sample configuration flexibility. This includes more readily achievable stability constraints for spectroscopies ( $\mu$ EXAFS), accommodating large sample sizes (up to meters), and provisions for customized environmental chambers.



## APPLICATIONS





- XFM will provide the NSLS-II user community an optimized beamline for studying the genetic control of metal ion uptake, transport and storage in plants relevant to agriculture and bioenergy.
- The only beamline in the world designed to directly support plant biochemistry, XFM will provide high-throughput and high-resolution whole-plant fCMT.
- fCMT provides a 3D non-invasive, spatially resolved and multi-elemental analysis technique that images the metal concentration of specific cell layers and organelles in plants as close to their natural state as possible.

- XFM's unique optical design, which allows it to maintain a small spot size at long working distances with high flux and at high energies, will provide an ideal platform for microfocused analysis of samples within environmental cells.
- For example, we are designing cells for XFM that allow users to analyze materials under  $\text{scCO}_2$  confinement ( $T > 31^\circ\text{C}$ ,  $P > 7.4\text{ MPa}$ ), while controlling  $\Delta P$ .
- XFM will utilize in-situ  $\mu$ XRF,  $\mu$ EXAFS and  $\mu$ XRD to quantify hydraulic/transport properties of brine films confined by  $\text{CO}_2$  under geologically relevant conditions.

- Imaging large samples at X-ray microprobes is impractical, but XFM with its variable focus and collimation and advanced ultrafast, large solid-angle EDS detectors being developed at BNL, is ideal for these studies.
- XFM will allow for signals from 100's of elements to be read at milliseconds per pixel for bidirectional scans covering meters!
- Large format translation stages and environmentally-controlled hutches will allow for analysis of whole objects including panel paintings, sculpture, paleontological and archaeological materials.
- XFM will be the focus of such work for museum scientists in the northeastern U.S.

## WORLD-LEADING MICROFOCUSED EXAFS SPECTROSCOPY

Fig. 1. Optical microscope image of *Thlaspi* seed cryosection,  $\mu$ XRF images of Ni and S in seed epidermal and palisade cells, and Ni K-edge  $\mu$ EXAFS spectra from single epidermal cell.



- XFM's unique optical design is optimized to enhance NSLS-II's source stability to provide unmatched  $\mu$ EXAFS data quality. Specifications include:
  - Use of a toroidal focusing mirror to illuminate a virtual source on the experimental floor that will be collimated, monochromated and imaged with KB microfocusing mirrors.
  - $\mu(E)$  oscillations measured to  $<0.1\%$  of the edge-step signal for transition metals at  $>10\text{ ppm}$  concentration in a  $1\text{-}10\text{ }\mu\text{m}$  spot.
  - Focused beam stability of  $<5\%$  of the beam area while scanning a  $1\text{ keV}$  energy range.
  - Fixed-exit DCM with user-selectable Si(111) and Si(311) crystal-pairs on a broadband source, will permit general user's to seamlessly and efficiently configure the beamline for spectroscopy in the  $4\text{-}20\text{ keV}$  range.