

***LABORATORY CHARACTERIZATION OF A SIZE-RESOLVED CPC
BATTERY TO INFER THE COMPOSITION OF FRESHLY FORMED
ATMOSPHERIC NUCLEI***

Kuang, C., Kangasluoma, J., Wimmer, D., Lehtipalo, K., Wang, J.,
Kulmala, M., and Petäjä, T.

Extended Abstract

For presentation at the
19th International Conference on
Nucleation and Atmospheric Aerosols
Fort Collins, CO
June 24-28, 2013

Atmospheric Sciences Division/Environmental Sciences Dept.

Brookhaven National Laboratory

**U.S. Department of Energy
Office of Science**

Managed by
Brookhaven Science Associates, LLC
for the United States Department of Energy under
Contract No. DE-AC02-98CH10886

DISCLAIMER

This report was prepared as an account of work sponsored by an agency of the United States Government. Neither the United States Government nor any agency thereof, nor any of their employees, nor any of their contractors, subcontractors, or their employees, makes any warranty, express or implied, or assumes any legal liability or responsibility for the accuracy, completeness, or any third party's use or the results of such use of any information, apparatus, product, or process disclosed, or represents that its use would not infringe privately owned rights. Reference herein to any specific commercial product, process, or service by trade name, trademark, manufacturer, or otherwise, does not necessarily constitute or imply its endorsement, recommendation, or favoring by the United States Government or any agency thereof or its contractors or subcontractors. The views and opinions of authors expressed herein do not necessarily state or reflect those of the United States Government or any agency thereof.

Laboratory Characterization of a Size-Resolved CPC Battery to Infer the Composition of Freshly Formed Atmospheric Nuclei

Chongai Kuang^a, Juha Kangasluoma^b, Daniela Wimmer^{b,c}, Katrianne Lehtipalo^{b,d}, Jian Wang^a, Markku Kulmala^b, Tuukka Petäjä^b

^a *Atmospheric Sciences Division, Brookhaven National Laboratory, Upton, New York, 11789, USA*

^b *Department of Physics, P.O. Box 64, 00014, University of Helsinki, Helsinki, Finland*

^c *Institute for Atmospheric and Environmental Sciences, Goethe University of Frankfurt, Frankfurt am Main, Germany*

^d *Airmodus, Ltd., Finland, Gustaf Hällströmin katu 2 A, 00560 Helsinki, Finland*

Abstract. A size-resolved condensation particle counter battery (SR-CPCb) was developed to infer the composition of freshly nucleated aerosol down to 1 nm in diameter. This was accomplished by exploiting the strong dependence of CPC counting efficiency on particle composition and CPC working fluid. The SR-CPCb consists of CPCs that are optimized for sub 2 nm particle detection using the following working fluids: diethylene glycol, water, and butanol. Laboratory characterizations of the SR-CPCb were performed for challenge aerosols of various compositions down to 1 nm in diameter.

Keywords: CPC calibration, heterogeneous nucleation, particle activation.

PACS: 92.60.Mt, 82.33.Tb

INTRODUCTION

Atmospheric particle nucleation is an important environmental nano-scale process, with climate models indicating that nearly half of the global cloud condensation nuclei (CCN) may be formed from freshly nucleated particles. However, our understanding of atmospheric nucleation and its influence on climate is limited since few direct measurements have been made of either the nucleation rate or the chemical composition of the freshly formed clusters, both of which are necessary to constrain the mechanism and develop a nucleation model. In our study, we have developed an instrument to infer the size-resolved composition of freshly formed atmospheric particles, addressing a key knowledge gap in the composition of nucleated aerosol below 3 nm. This size-resolved condensation particle counter battery (SR-CPCb), is an extension of the CPCb concept developed and used to infer the composition of nucleation mode particles between 2 and 9 nm [1, 2]. In the SR-CPCb, the CPCb measurement principle has been extended down to 1 nm by combining a conventional nanoparticle mobility classifier with a CPCb composed of CPCs designed for the detection of sub 2 nm aerosol [3-7]. By sampling mobility-classified particles, the size resolved CPCb (SR-CPCb) eliminates the strong dependence of CPC response on particle size and charge below 2 nm [4-10]. Therefore, any measured differences in

CPC response are then attributed solely to composition-specific interactions between the particle and the various working fluids.

EXPERIMENTAL

The SR-CPCb is composed of a nanoparticle mobility spectrometer that has been optimized for the sampling and mobility classification of sub 3 nm particles [9]. The particle detector is composed of a CPCb containing the following CPCs with their accompanying working fluids: a Particle Size Magnifier [3] (diethylene glycol), a DEG-UCPC [5] (diethylene glycol), a TSI 3786 UCPC (water), a TSI 3025A UCPC [5] (butanol), and a pulse height UCPC [6] (butanol). The experimental schematic for determining the size, charge, and composition dependent particle detection efficiencies of the CPCb is described in [10]. The following methods were used to generate challenge aerosols for the characterization of the CPCb: (1) electrospray generation of molecular ion mobility standards, (2) evaporation of solid sodium chloride via tube furnace, (3) evaporation of ammonium sulfate via tube furnace, (4) tungsten oxide formation via wire generation, and (5) candle aerosol generation. Presented results will include the particle composition dependent response of the SR-CPCb and size-dependent transmission efficiencies associated with the nanoparticle mobility spectrometer.

ACKNOWLEDGMENTS

This work was supported by the US Department of Energy (DOE) Atmospheric System Research Program (Office of Science, OBER) under contract DE-AC02-98CH10886, US DOE Grant Number DE-FG-02-05ER63997, the Academy of Finland (project nos. 1251427, 1139656, 141451, 133872), the Finnish Center of Excellence (project no. 1118615), and the European Research Council (ATMNUCLE).

REFERENCES

1. Riipinen, I., et al., *Applying the Condensation Particle Counter Battery (CPCB) to study the water-affinity of freshly-formed 2–9 nm particles in boreal forest*. *Atmos. Chem. Phys.*, 2009. **9**: p. 3317-3330.
2. Kulmala, M., et al., *The condensation particle counter battery (CPCB): A new tool to investigate the activation properties of nanoparticles*. *Journal of aerosol science*, 2007. **38**(3): p. 289-304.
3. Vanhanen, J., et al., *Particle Size Magnifier for nano-CN detection*. *Aerosol Science and Technology*, 2011. **45**(4): p. 533-542.
4. Iida, K., M.R. Stolzenburg, and P.H. McMurry, *Effect of working fluid on sub-2 nm particle detection with a laminar flow ultrafine condensation particle counter*. *Aerosol Science and Technology*, 2009. **43**(1): p. 81-96.
5. Kuang, C., et al. *Modification of Laminar Flow Ultrafine Condensation Particle Counters for the Enhanced Detection of 1 nm Condensation Nuclei*. *Aerosol Science and Technology*, 2012. **46**, 309-315 DOI: 10.1080/02786826.2011.626815.
6. Sipilä, M., et al., *Laboratory Verification of PH-CPC's Ability to Monitor Atmospheric Sub-3 nm Clusters*. *Aerosol Science and Technology*, 2009. **43**(2): p. 126-135.
7. Wimmer, D., et al., *Performance of diethylene glycol based particle counters in the sub 3 nm size range*. *Atmospheric Measurement Techniques Discussions*, 2013.

8. Winkler, P.M., et al., *Heterogeneous nucleation experiments bridging the scale from molecular ion clusters to nanoparticles*. *Science*, 2008. **319**(5868): p. 1374.
9. Jiang, J., et al., *Electrical Mobility Spectrometer Using a Diethylene Glycol Condensation Particle Counter for Measurement of Aerosol Size Distributions Down to 1 nm*. *Aerosol Science and Technology*, 2011. **45**(4): p. 510 - 521.
10. Kangasluoma, J., et al., *Remarks on Ion Generation for CPC Detection Efficiency Studies in Sub 3 nm Size Range*. *Aerosol Science and Technology*, 2013: p. null-null.