

Stimulated X-ray Emission Spectroscopy with hard X-rays

Thomas Kroll

SLAC National
Accelerator Laboratory

Acknowledgement

SLAC

SLAC

Uwe Bergmann
Clemens Weninger
Y. Zhang
Roberto Alonso-Mori
Agostino Marinelli
Alberto Lutman
Marc W. Guetg

CXI Instrument at LCLS

Sebastien Boutet
Andy Aquila
Jake Koralek
Dan DePonte

Max Planck, Hamburg, Germany

Nina Rohringer
Laurent Mercadier
Vinay P. Majety

Lawrence Berkeley National Lab

Vittal Yachandra
Junko Yano
Jan Kern
Franklin Fuller





MAX-PLANCK-GESELLSCHAFT



Funding: DOE and NIH

X-Ray Emission Spectroscopy

Level Diagram



Spin/oxidation state of transition metals


Valence orbitals: ligand type, structure, covalency, ligand protonation, etc

Bergmann et al, *J Synchr Rad.* 8, 199 (2001)

Stimulated X-ray Emission Spectroscopy

Stimulated (non-linear) emission:

Spontaneous (conventional) emission




- Low photon flux density (only a few photons at a time in the sample)
- Only spontaneous emission in 4π direction

Stimulated X-ray Emission Spectroscopy

Stimulated (non-linear) emission:


Amplified spontaneous emission (ASE)

Pump:
(> abs. energy)



- Very high flux density ($\sim 10^{11}$ ph / 10 fs / 100x150 nm²)
- Creation of population inversion
- Non-linear effects become dominant
- Cascade like decay


Experimental Setup



- Highly focused 10 fs beam ($100 \times 150 \text{ nm}^2$)
- Flat crystal analyzer in forward direction



- Can stimulated XES be applied to transition metal solutions?
- What is fundamentally new information provided by stimulated XES techniques that is not achievable with other techniques?
- Can stimulated XES be used for chemically sensitive X-ray emission lines? Is the chemical information preserved?
- Can we stimulate the weaker $K\beta$ emission?

Single shot spectra ($K\alpha$ -ASE)





SLAC

Single shots:



- Clean curve in single shot
- Only the $K\alpha_1$ visible
- FWHM < normal $K\alpha$ width



Lasing in MnCl_2 solution



Proof of Lasing:


- Exponential increase
- Linear gain regime below ~ 2 mJ
- Saturation reached > 2 mJ

Peak Position and Width





- Low to mid high photon numbers:
 - Constant broadening and position
- High photon numbers:
 - Spectral broadening
 - Shift to lower energies
- Variations through beam position, temporal shape, lasing condition



Peak Position and Width



- Line broadening due to saturation effects
- Additional final states become visible



Peak Position and Width



- Mn $K\alpha_1$ life-time broadening: 1.48 eV (Krause and Oliver, 1979)
- Lowest observed S-XES peak width: < 1.0 eV
- Darwin width of Si (111): 0.77 eV
- Lowest S-XES lifetime broadening: < 0.5 eV


Stimulated X-ray Emission Spectroscopy

Stimulated (non-linear) emission:

Amplified spontaneous emission (ASE) + seeded stimulated emission


Pump:
> abs. energy

Seed:
 $K\beta$ energy




- $K\beta$ Seed pulse outruns $K\alpha$ emission → seeded stimulated $K\beta$ emission
- Requirement: Overlap in time, space and energy

Experimental Setup



- Highly focused 10 fs beam ($100 \times 150 \text{ nm}^2$)
- Use part of the undulators for seed pulse
- Two flat crystals in series in forward direction for simultaneous detection of $K\alpha$ and $K\beta$

Seeded Stimulated K β emission



Pump:
> abs. energy

Seed:
K β energy

Biggest issue:

- Seed pulse and stimulated emission are at the same energy
- How to separate the stimulated signal from the SASE seed pulse?

Conclusions

- Strongly stimulated XES in Mn solutions has been observed
- Gain curves for X-ray lasing have been measured
- Gain narrowing to less than the natural line width
- Chemical information appears to be preserved!
- Seeded stimulated $K\beta$ emission observed

Future



- Improved diagnostics:
 - especially in time domain
 - Seeded beam diagnostic

Stimulated X-ray Emission has been observed

Past achievements:


Soft X-rays

- Stimulated soft X-ray Ne laser
(Rohringer et al, *Nature* **481**, 488 (2012), LCLS)
 - Ne gas
- Stimulated X-ray Emission in Silicon L Lines (70-100 eV)
(Beye et al, *Nature* **501**, 191 (2013), LCLS)
 - Si





Hard X-rays

- Stimulated hard X-ray Fe laser
(Yoneda et al., *Nature* **524**, 446 (2015), SACLA)
 - Fe foil



High number of photons:

Varying pulse energy:




- Increase of K α ₂ for high fluences
- Can be explained by final state effects

Chemical Sensitivity

Compare solution spectra:

- Stimulated data retain the expected chemical shift
- Need to address stimulated photon / matter interaction in more detail



Experimental Setup

X-ray Free Electron Laser: LCLS



LCLS @ SLAC Linear Accelerator Laboratory

- > 3 km long building
- Electrons are accelerated
- Undulators create photons
- Pulsed and highly focused beam
- Large number of photons per pulse:
 - 10^{11} per pulse of 10 fs