



Resonant Inelastic X-ray Scattering at SIX of NSLS II: An update on the beamline status and its performances

Valentina Bisogni

Brookhaven National Laboratory

Acknowledgements:



SIX Beamline Team



Andi Barbour
Ignace Jarrige
Jonathan Pellicari

Optics Scientist
Joe Dvorak

Soft X-ray Scattering & Spectroscopy Group



Claudio Mazzoli
Daniel Mazzone
Andrew Walter
Stuart Wilkins
....

Mechanical Engineers

Stephen Antonelli
Daniel Bacescu
Lukas Lienhard
Bill Leonhardt
Yi Zhu

Control Engineers

John Sinsheimer
Jun Ma
Wayne Lewis

Technicians

Larry Fareira
Andy Mingino
Joe Sullivan

And many more ...


Outline:

A. Soft Inelastic X-ray Scattering Beamline SIX at NSLS -II

- Beamline overview
- Status update: flux, BL energy resolution, beam size at the sample, spectrometer resolution

B. Collective excitations in model devices

- Scope of the research project
- First data measured at SIX



2-ID

SIX



Design Parameters for Resonant Inelastic X-ray Scattering at SIX

Energy: soft x-rays, 200 – 2000 eV

Polarization of incoming light:
LV, LH, Circular


Operation modes:
Medium Resolution ($R \sim 15000$)
High Resolution ($R \sim 35000$)
Ultra High Resolution ($R \sim 70000$)

Scattering angle:
Continuous in vacuum rotation
from 38 deg to 150 deg



Beamline open to user proposals:
Next Call: September 30th

<https://www.bnl.gov/ps/beamlines/>


SIX Beamline overview* : 105m of Photon Delivery System



Centurion Spectrometer: M5 collecting mirrors




Two CCD sensors imaging each one of two beams from the parabolic mirror



- M5 alignment requires a μm -sensitive alignment procedure
- It is fundamental to have the M5 focal point, the center of rotation of the 2 θ optics wheel and the M4 focal spot all at the same point!


Review of Scientific Instruments **87**, 115109 (2016)

Centurion Spectrometer: M6/M7/GR focalize and analyze energy




- Hettrick-Underwood scheme + Dvorak's mirror (M7) keep outgoing beam angle fixed at -1° at all energies
- CCD height is fixed

M6(M7) pitch scan at fixed GR pitch





GR(M7) pitch scan at fixed M6 pitch



Review of Scientific Instruments **87**, 115109 (2016)

SIX Beamline: Flux




- GR 500 lmm^{-1} provides 6 x more flux than GR 1800 lmm^{-1}
- Flux linear with exit slit vertical gap: **medium res mode** (GR 500 lmm^{-1} , 30 μm es vg) **~40x** more flux than **ultra high res mode** (GR 1800 lmm^{-1} , 5 μm es vg)

SIX Beamline: Energy Resolution for GR 500 Imm⁻¹

- XAS transmission measurement of reference gases

GR 500 Imm⁻¹ – Resolving Power ~ 20000



- ✓ Calibrate the PGM energy
- ✓ Estimate the energy resolution
- ✓ Identify the optimum C_{ff} value

BL Energy: 244 eV
 Measured Res: 9 meV
 Simulated Res: 7.2 meV

400 eV
 22.4 meV
 13.9 meV


867 eV
 44.9 meV
 42.4 meV

=> R ~ 19000 – 27000

SIX Beamline: Energy Resolution for GR 1800Imm⁻¹

- XAS transmission measurement of reference gases

GR 1800 Imm⁻¹ – Resolving Power ~ 100000




BL Energy: **400 eV**
 Measured Res: **0 meV**
 Simulated Res: **4 meV**

867 eV
0 meV
8.7 meV

- ✓ Calibrate the PGM energy
- ✗ Estimate the energy resolution
- ✗ Identify the optimum C_{ff} value

- PGM is calibrated when the same energy goes through the exit slit during a C_{ff} scan




- Fine tuning of C_{ff} optimal value is done using the spectrometer


SIX Beamline: Beam Size at the Sample

- New M4 refocusing mirror from JTEC installed in April 2019


EXIT SLIT VG 5 μm




EXIT SLIT VG 10 μm




Vertical Beam Profile vs Exit Slit VG



Vertical Focus Size





Performance close to ideal while using a 2 mm vertical gap on the baffle slits prior to the mirror



Centurion Spectrometer: Current Performances



BL GR 500 Imm^{-1} / SP GR 2500 Imm^{-1}

Resolving Power ~ 21000




BL GR 1800 Imm^{-1} / SP GR 2500 Imm^{-1}



Resolving Power ~ 42000





Centurion Spectrometer: Scattering Arm Rotation



- Spectrometer arm motion system fully implemented and commissioned.
- Sample chamber is equipped with a discrete number of fixed exit ports.
- Long term solution: fully in vacuum rotation with Triple Rotating Flange



Triple Rotating Flange Concept:





Thanks for your attention!