

HAXPES on the beamline ID32 at the ESRF and its application to buried interfaces in exchange bias heterostructures



Blanka Detlefs, Lionel André, Helena Isern, Yanyu Mi, Lucien Petit, Michel Renier, Jérôme Roy, and Jörg Zegenhagen

European Synchrotron Radiation Facility, 6 rue Jules Horowitz, F-38043, Grenoble, France

Tien-Lin Lee

Diamond Light Source Ltd., Rutherford Appleton Laboratory, Chilton-Didcot OX11-0QX, UK

Sara Laureti and Dino Fiorani


C.N.R. - ISM Via Salaria, Monterotondo Stazione, I-00016, Roma, Italy

Science at the beamline ID32

- Surface X-ray diffraction: determination of surface structure, electronic properties and solid/electrolyte interfaces. In-situ electrochemistry
- X-ray standing waves (combined with XPS, SXRD): site specific XPS, imaging of surface structures and adsorbates
- Surface spectroscopy - SEXAFS; high energy resolution available for resonant scattering and absorption spectroscopy
- HAXPES

Beamline setup

- windowless beamline, 3 undulators (1 to ~35 keV)
- fixed-exit, liquid nitrogen cooled high-precision double crystal monochromator; with Si(111) crystals: photon energies 2.1 to 30 keV, energy resolution $\Delta E/E = 1.3 \cdot 10^{-4}$
- flux $> 10^{13} \text{ mm}^{-2}\text{s}^{-1}$
- 2nd set of multilayer crystals installed for possible experiments below 2 keV
- main mono 2nd crystal piezo control: for vertical beam position feedback (diffraction experiments)
- flat mirror (Ni, Pd and SiO₂ coating) for higher harmonics rejection
- parabolic Be compound refractive lenses (2:1 focusing) \rightarrow beam-size at sample $30 \times 30 \mu\text{m}^2$
- Si Fresnel zone plates (30:1 focusing) \rightarrow beam-size at sample $4 \times 15 \mu\text{m}^2$ at certain energies




High energy resolution: post-monochromator

- high energy resolution at variable energy
- two axes, presently with Si(111), Si(220) and Si(311) channel-cut crystals
- dispersive arrangement (+/-, +/-) for high energy resolution down to $\Delta E/E = 10^{-6}$
- backscattering condition (fixed energies) to match photon beam divergence
- asymmetrically cut crystals for flexible energies (in preparation)


HAXPES

- Phoibos225 electron analyzer by SPECS:
 - for electron energies up to 15 keV
 - 2D delay-line detector (or CCD): excellent signal-to-noise ratio
 - different modes of operation
 - UPS, XPS and HAXPES
 - angular or spatial resolution in the non-dispersive direction
 - scanning or imaging mode
 - ultra high energy resolution (specifications: in UPS $< 1 \text{ meV}$, in HAXPES $< 15 \text{ meV}$)
- HAXPES chamber to be equipped with
 - 4-circle kappa goniometer
 - preparation chamber (early 2010)
 - He flow cryostat (summer 2009)




HAXPES commissioning (May 2008 – May 2009)

- HV power supply operational up to 15 kV
- data acquisition up to 13 kV successful




- energy resolution tests:
 - at $E = 8.0 \text{ keV}$, using Si(444) post-monochromator (one axis only): $\Delta E_{\text{phot}} \sim 40 \text{ meV}$, at room temperature
 - $\Delta E_{\text{tot}} = 102(6) \text{ meV}$
 - $\rightarrow \Delta E_{\text{analyzer}} < 30 \text{ meV}$




Exchange bias

- exchange coupling at the interface between a ferromagnetic (FM) and an antiferromagnetic (AFM) material
- shift of the hysteresis loop (exchange bias field, HE) is observed when the sample is cooled in a magnetic field through the Néel temperature: induced unidirectional anisotropy that arises from the direct exchange coupling at the FM/AFM interface
- the exchange coupling usually leads also to an increase in coercivity (Hc), as found in both nanoparticles and bilayer EB systems
- wide technological applications (spin-valve systems, GMR read heads, tunneling devices)
- complexity of the interface microstructure
- factors: interface roughness, thickness of both FM and AFM phases, interface stress etc.
- control of growth mechanism
- detailed characterization of the interfacial region




Exchange bias bilayers

- fcc-CoPt(FM)/NiO (AFM) fabricated by PLD on Si(100) substrate (+ native oxide)
- different thickness of NiO (10- 50 nm) + 5 nm of CoPt: ablation of a composite Co and Pt target
- laboratory characterization: EDXD and SQUID measurement
- X-ray reflectivity measurements: example - sample with 11 nm NiO + 5 nm of CoPt: roughness (SiO₂/NiO) = 0.3 nm, roughness (NiO/CoPt) = 0.2 nm



Depth profiling using HAXPES

- probing depth increases with increasing photon energy: Ni 2p and Si 2s, 2p states discernible
- probing depth increases with increasing take-off angle
- Co 2p states: lower CoO signal (vs. metallic Co corresponding to Co bound to Pt) with increasing probing depth (both increasing E_{phot} , increasing take-off angle): CoO present mostly on the bilayer surface



Summary

- HAXPES operational at the beamline ID32 at the ESRF
- high energy resolution available at variable photon energy
- enlarged HAXPES possibilities in future - low temperature, X-ray standing wave technique, in-situ sample preparation and characterization, remote control operation, in-situ electrochemistry
- depth profiling of heterostructures of ~10 - 15 nm possible at energies above 5 keV: quantitative analysis difficult in heterostructures with (possibly) repeating chemical composition
- combination of HAXPES with reflectivity measurements desirable for valid structural model

Acknowledgment

We would like to thank SPECS and Surface Concept staff (in particular S. Maehl, A. Oelsner, C. Amicabile, S. Egner) for a continuous support in the commissioning of SPECS 225 electron analyzer.